

## MinION Mk1B (MN19362) Final report

16 Jun 23, 17:24 — 16 Jun 23, 17:49 · RU\_test · no\_sample · MN19362  
 Protocol run ID: cf0b707b-477e-454b-82c2-f532dae16fc5




[Run summary](#) | [Run configuration](#) | [Sequence output](#) | [Run health](#) | [Run log](#)

### Run summary

#### DATA OUTPUT

Estimated bases	Data produced
77.41 Mb	679.25 MB
Reads generated	Estimated N50
4.44 k	38.23 kb

#### RUN DURATION



[View unit abbreviations used in this report](#)

#### BASECALLING

Reads called	Bases called (min Q score: 8)	
90.74%	38.37 Mb	5.16 Mb
	Pass	Fail

### Run configuration

#### RUN SETUP

Flow cell type	FLO-MIN106
Flow cell ID	sim
Kit type	SQK-LSK109

#### RUN SETTINGS

Run limit	72h
Active channel selection	On
Pore scan freq.	1.5 hrs
Reserved pores	On
Minimum read length	200 bp
Read splitting	Off
Basecalling	Fast model, 450 bps
Modified basecalling	Off

#### DATA OUTPUT SETTINGS

FAST5 output	vbz_compress
FAST5 reads per file	4000
FASTQ output	gzip_compress
FASTQ reads per file	4000
BAM output	Off
Bulk file output	Off
Data location	/var/lib/minknow/data/.RU_test/no_sample/20230616_1324_MN19362_sim_cf0b707b

#### SOFTWARE VERSIONS

MinKNOW	23.04.5
Bream	7.5.9
Configuration	5.5.13
Guppy	6.5.7
MinKNOW Core	5.5.3

## Sequence output

### READ LENGTHS · OUTLIERS REMOVED


The read length graph shows the total number of bases vs the read length. The longest 1% of strands are classified as outliers, and excluded to allow focus on the main body of data.

Legend

■ Basecalled — Estimated

Estimated N50  
38.23 kb

% Basecalled  
90.74%



### OUTLIERS

The longest 1% of strands are classified as outliers, and aggregated into groups to show their relative amounts.

Read length (kb)	Aggregated reads (kb)
392 - 396	383.7


### CUMULATIVE OUTPUT

The cumulative output shows the total amount of bases or reads sequenced over time by your device.

#### Bases

Legend


— Estimated: Predicted total number of bases, prior to basecalling  
 — Passed: Bases equal to or above the quality score threshold.  
 — Failed: Bases below the quality score threshold.



#### Reads

Legend

— Total: Total number of reads, including passed, failed and skipped.  
 — Passed: Reads equal to or above the quality score threshold.  
 — Failed: Reads below the quality score threshold.  
 — Skipped: Reads that will not be basecalled. Post run basecalling is possible.




### QUALITY SCORE

The quality score is calculated as basecalling is performed on your device. Reads that fall below the minimum value of 8 will be classified as failed reads. You can alter the accepted minimum quality score in MinKNOW.

Legend

- Mode  
The most frequent quality score of reads in the run.
- Spread  
The spread of quality scores, found by calculating full width half maximum.
- Min. quality score  
Minimum quality score to be accepted as a passed read.



Troubleshooting

Quality score low

This can be due to the translocation speed being out of the accepted range, which can correlate to low quality scores. If you see that the translocation speed is out of the accepted range in the below graph, please see the Flow Cell refuelling page linked [here](#) for further troubleshooting.


### Run health

#### PORE ACTIVITY

The Pore activity graph shows the performance of your sample as it is being sequenced during a run.

Legend

- Sequencing  
Pore currently sequencing
- Adapter  
Pore currently sequencing adapter
- Pore available  
Pore available for sequencing
- Unavailable  
Pore unavailable for sequencing
- Active feedback  
Channel ejecting analyte
- No pore  
No pore detected in channel
- Out of range-high  
Current is positive but unavailable for sequencing
- Out of range-low  
Current is negative but unavailable for sequencing
- Multiple  
Multiple pores detected. Unavailable for sequencing.
- Saturated  
The channel has switched off as current levels exceed hardware limitations
- Zero  
Pore currently unavailable for sequencing
- Channel disabled  
Channel is disabled and awaiting another pore scan
- Unclassified  
Pore status unknown



Troubleshooting

General


Some commonly seen issues are excess pores classified as Recovering, Open Pore, or Free Adapter. To find out what advice is applicable for your run, visit the [user guide](#).

### PORE SCAN

A Pore scan is performed at configurable time intervals to determine the current status of pores within channels on a Flow Cell. For this run a Pore scan is performed every 1.5 hrs.

Legend

- **Pore available**  
Pore in channel available for sequencing
- **Reserved pore**  
Pore in reserve, will return to available when required
- **Unavailable**  
Pore inhibited from sequencing
- **Saturated**  
Possible contamination in the sample
- **Zero**  
No current is passing through this pore, possibly due to bubbles on the membrane
- **Inactive**  
Pore no longer suitable for further sequencing



#### Troubleshooting

**High proportion Unavailable**  
Possible contaminants in library blocking the pore. Consider using the Flow Cell Wash Kit, and reloading a library.


**High proportion Inactive**  
If localised to one area of the Flow Cell, this could indicate that an air bubble has been introduced during the flushing/loading steps. If inactivity is spread across the Flow Cell this could be caused by improper loading of the library, please refer to the [user guide](#) for further support.

### TRANSLOCATION SPEED

The translocation speed is the rate at which DNA/RNA travels through pores as it is being sequenced.

Legend

- Median
- - - 75% quartile
- - - 25% quartile
- Accepted range



#### Troubleshooting

**Low speed**  
Check that the Flow Cell is within the target temperature range.


**Note**  
Low-quality and short reads are not included in this graph.

### TEMPERATURE

The temperature of the Flow Cell over the run time.

Legend

- Measured
- Target



#### Troubleshooting

**Out of range**  
Check that the Flow Cell is correctly seated and firmly pushed down into the device. Ensure ambient temperature is always within the specified range for your device in the [user guide](#).

Air flow should be good but not excessive. Excessive amounts of cool air blowing on the device could prevent it from reaching target temperature.

## Run log

### SYSTEM MESSAGES

System messages are a record of the events that occurred in the time covered by this report.

Please view the original run report, which would have been generated at the end of the run, to see the system messages.

---

### UNIT ABBREVIATIONS

Byte	B
Kilobyte	KB
Megabyte	MB
Gigabyte	GB
Terabyte	TB

Base	b
Kilobase	kb
Megabase	Mb
Gigabase	Gb
Terabase	Tb

Minutes	mins
Hours	hrs

Generated using Run Report Template v.5.5.6